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RevitalVision

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V1

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V1.0.0

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The login page features a blue header with the text "针对大脑的训练，使您看得更清楚". Below the header is a "登录" (Login) button with a user icon. Underneath is a form with two input fields: the first contains the username "tt003" and the second contains masked characters "*****". A "登录" (Login) button is positioned below the password field. At the bottom, there is a "新用户?" (New User?) link and a "尝试免费演示版本" (Try Free Demo Version) button.

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The system homepage has a blue header with the "REVITALVISION" logo and navigation links: "训练" (Training), "教学" (Teaching), "使用说明" (Usage Instructions), "关于" (About), and "联系" (Contact). A language dropdown menu is set to "Chinese". On the left, a sidebar shows the user "tt003" and navigation options: "我的训练" (My Training), "我的进展" (My Progress), and "校准" (Calibration). The main content area displays a welcome message "欢迎您, tt003,".

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The calibration interface is divided into two sections: "屏幕尺寸校准" (Screen Size Calibration) and "亮度校准" (Brightness Calibration). The "屏幕尺寸校准" section includes instructions to adjust the screen width to 10 cm and height to 5 cm using a ruler. It features a blue bar for width adjustment and a red bar for height adjustment, both with a red arrow icon. Below these are "完成" (Complete) and "取消" (Cancel) buttons. The "亮度校准" section shows a grayscale calibration chart with two rows of numbers: "3 6 0 2 9 0" and "5 0 4 5 9 9". Below the chart are two input fields for the first and second rows of numbers, and a "确定" (Confirm) button.

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